KERN **EMS3000-2**

Entry level laboratory balance with large weighing plate





| Measuring system | |
|---|----------------------|
| Adjustment options: | External calibration |
| Eccentric loading at 1/3 [Max]: | 50 mg |
| Linearity: | 50 mg |
| Readability [d]: | 10 mg |
| Recommended adjusting weight: | 3 kg (F2) |
| Repeatability: | 20 mg |
| Resolution: | 300000 |
| Stabilisation time under laboratory conditions: | 3 s |
| Tare range: | 3 kg |
| Warm up time: | 2 h |
| Weighing capacity [Max]: | 3 kg |
| Weighing system: | Strain gauge |
| Weighing units: | g dwt oz ozt |
| Display | |
| Display digit height (large): | 2,500 cm |
| Counting | |
| Counting resolution: | 150000 |
| Mininum piece weight at piece counting (Laboratory: | 20 mg |
| Mininum piece weight at piece counting (Normal): | 200 mg |
| Functions | |
| Counting function: | yes |
| Percentage determination: | yes |
| Percentage determination: | yes |

| Environmental conditions | |
|---------------------------------|-------------------------|
| Maximum humidity: | 80 % |
| Maximum operating temperature: | 35 °C |
| Minimum ambient temperature: | 5 °C |
| Minimum humidity: | 40 % |
| Power supply | |
| Battery: | 9 V Block |
| Input voltage: | 220 V - 240 V AC 50 Hz |
| Mains adapter: | Mains adapter external |
| Mains adapter/adapter included: | UK EURO CH US |
| Operating time: | 40 h |
| Rechargeable battery: | Rchrg. battery optional |
| Services (optional) | |
| DAkkS Certificate: | 963-127 |
| Category | |
| Brand: | KERN |
| Category: | Balances |
| Product Group: | School balance |
| Packaging & shipping | |
| Delivery: | 24 h |
| Dimensions packaging (WxDxH): | 310 x 234 x 164 mm |
| Gross weight: | 1,980 kg |
| Net weight: | 1,575 kg |
| Shipping method: | Parcel service |
| Construction | |
| Casing material: | Plastic |
| Dimensions housing (WxDxH): | 200 x 280 x 48 mm |
| Material weighing plate: | Plastic |
| Weighing surface (WxD): | 175 x 190 mm |
| Piktogramme | |
| STANDARD | |
| | |



OPTION

DAkks +3 days